



# i-STAT Alinity v Gamme des tests

			Bilans complets	Gaz du sang		Electrolytes		Tests unitaires		
			CHEM8+	EC8+	CG4+	CG8+	E3+	6+	G	Crea
Hématologie	Hématocrite	Hct	●	●		●	●	●		
	Hémoglobine*	Hgb	●	●		●	●	●		
Biochimie	Taux d'azote uréique	BUN	●	●				●		
	Créatinine	CRE	●							●
	Calcium ionisé	iCa <sup>++</sup>	●			●				
	Glucose	GLU	●	●		●		●	●	
Electrolytes	Chlorure	Cl <sup>-</sup>	●	●			●	●		
	Sodium	Na <sup>+</sup>	●	●		●	●	●		
	Potassium	K <sup>+</sup>	●	●		●	●	●		
Acido-basiques	Potentiel Hydrogène	pH		●	●	●				
	Pression partielle de dioxyde de carbone	pCO <sub>2</sub>		●	●	●				
	Bicarbonate*	HCO <sub>3</sub>		●	●	●				
	Dioxyde de carbone total*	tCO <sub>2</sub>	●	●	●	●				
	Trou anionique*	AG	●	●						
	Excès de base*	BE		●	●	●				
Gaz du sang	Pression partielle d'oxygène	pO <sub>2</sub>			●	●				
	Saturation en oxygène*	so <sub>2</sub>			●	●				
Spécialité	Lactate	LAC			●					

\*Valeurs calculées